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
- |   |                                  |
|---|----------------------------------|
| 1 Reclining surface frame                               | 10 Hand switch                   |
| 2 Lower leg rest  | 11 Fastening latches (4)         |
| 3 Motor drive, foot section                             | 12 Base frame                    |
| 4 Thigh rest  | 13 Pedestals (4)                 |
| 5 Motor drive, head section<br>(not visible in picture) | 14 Connecting socket             |
| 6 Backrest  | 15 Transformer unit              |
| 7 Mattress brackets (4)                                 | 16 Motor drive height adjustment |
| 8 Trapeze pole sockets (2)                              | 17 Wiring harness                |
| 9 Frame connectors (2)                                  | 18 Guide bars (2)                |

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**Please notice:**

In the text of these instructions square brackets and bold cypers = **[00]** refer to the parts of the slatted frame.

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Burmeier GmbH & Co. KG

## PREFACE

Dear customer,  
Thank you for the confidence you have shown in BURMEIER.

The vertically adjustable slatted frame LIPPE IV is factory-checked for electrical safety and function and has left our factory in faultless condition.

Please take time to read this manual carefully. This applies especially to the assembly phase but also to daily use.

It is also a practical reference guide. Keep the instructions in a handy place.

These instructions provide information for the operator and user regarding the convenient handling and safe operation.

We wish you and your users much success in care-giving and we are convinced that our product will make a positive contribution.

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*Typ 178629*

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## GENERAL NOTES

The slatted frame LIPPE IV is available in various models.

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#### **Before putting the slatted frame into operation for the first time:**

■ Read these instructions carefully in order to avoid damage or incorrect operation.

■ Clean and disinfect the slatted frame before first use.

The slatted frame LIPPE IV complies with all requirements of the guidelines 93/42/EWG for medical products. It is classified as a dynamic medical product of the class I in accordance with "Medizinproduktegesetz (*German abbreviation: MPG § 13, Medical Products Act*)."

Observe your obligations as operator in accordance with medical devices "Medizinprodukte-Betreiberverordnung (*German abbreviation: MPBetreibV, Operators of Medical Products Ordinance*)" in order to ensure a dependable and safe operation of this medical device without endangering patients, users and third parties.

The slatted frame LIPPE IV is checked by an independent examining institute. As with all technical, electrical devices improper handling can lead to damage and/or injury.

**This manual contains safety guidelines, which must be observed. All persons who work with the slatted frame LIPPE IV must**

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familiarize themselves with these instructions and follow the safety guidelines.

## DEFINITIONS OF INVOLVED PERSONS

In this manual the following persons are involved:

### **Operator**

Operators (e.g.: Medical Supply Houses, Dealers and representatives, health insurance companies) are natural persons or legal entities, who utilize the slatted frame LIPPE IV or authorize its use. The operator is responsible for instructing the user in the proper use of the unit.

### **User**

Users are persons who as a result of training and experience are authorized to assemble, perform maintenance and/or operate the slatted frame LIPPE IV. Furthermore, the user can recognize and avoid possible dangers, and judge the clinical condition of the patient.

### **Patient**

In this manual a patient is an infirm or handicapped person who reclines in a bed with this slatted frame.

### **Specialist staff**

Specialist staff are employees of the operator who are authorized as a result of their education and training to deliver, assemble, dismantle and transport the slatted frame. Furthermore, they are familiar with the regulations for cleaning and disinfecting the unit.

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## SAFETY INSTRUCTIONS

The slatted frame LIPPE IV represents the latest technological development at the time of delivery and is checked by an independent examining institute.

■ Do not use the slatted frame LIPPE IV if it is damaged.

### MEANING OF THE SAFETY SYMBOLS

In these instructions the following safety symbols are used:

#### Warning about injuries to persons

Dangerous voltage. Life Threatening.

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General danger. Injury or life hazard.

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#### Warning about property damages

Possibility of damage to drive, material or environment.

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#### Other symbols

Useful tip. For easier operation or better understanding of the unit.

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**The safety symbol used does not replace the text of the safety note. Therefore read the safety instructions and follow them exactly!**

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## SAFETY INSTRUCTIONS FOR OPERATORS

- Remember your responsibilities in accordance with the MPBetreibV in order to ensure the safe operation of this unit without endangering patients, users and third parties.
- Before first use, instruct every operator in the safe operation of the slatted frame in accordance with these safety instructions, which must be supplied along with the unit.
- Call attention to the dangers of improper use of the unit, especially in regards to the electrical drives and sideguards.
- The slatted frame should only be operated by properly trained personnel.
- Make sure that representatives are also familiar with the operation of the slatted frame.
- During long term use it is recommended that the slatted frame be inspected for function and visible damage at regular intervals (Recommendation: yearly) (see chapter 6).
- Connect the slatted frame directly to an electrical outlet. Do not use extension cords and/or multiple outlet extension cords.
- When attaching other equipment (e.g. compressors for positioning systems) make sure that they are mounted and can function in a secure and safe manner.
- Do not put multiple outlet electrical extensions under the slatted frame. Leaking liquids can cause a fire hazard.
- Pay special attention to:  
Secure locating of all wiring, cables, tubes, etc.  
▶ See page 12: *Intended use*

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▶ See page 26: *Safety note for the electric drive system*

For further information please contact the manufacturer of the accessories or BURMEIER directly.

■ Make sure that your personnel follow the safety instructions.

### SAFETY INSTRUCTIONS FOR USERS

■ Let the operator instruct you in the safe use of the slatted frame.

■ Before each use make sure that the slatted frame is in proper and faultless condition.

■ Take care that no obstacles, such as furniture or slanted ceilings interfere with the adjustment functions.

■ Pay attention that when using additional electrical components, such as patient lifts, reading lights or compressors for positioning systems that their electrical cords do not get entangled or damaged by the moving parts of the slatted frame.

■ Make sure that when attaching other appliances (e.g. compressors for positioning systems) that a safe mounting and function is guaranteed.

■ Do not put multiple electrical outlets under the slatted frame. Leaking liquids can be a fire hazard.

■ Pay special attention to securely locating all wiring, cables, tubes, etc.

For further information please contact the manufacturer of the accessories or BURMEIER directly.

■ Take the slatted frame out of operation if damage or a malfunction is suspected:

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- ▶ Immediately unplug the electrical plug from the wall outlet;
- ▶ Clearly mark the slatted frame as " DEFECTIVE ";
- ▶ Promptly notify the appropriate operator.

A checklist for assessing the proper condition is located in Chapter 6.

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- Position the power cable in such a way that during normal operation of the slatted frame the cable will not be stretched, driven over or interfere with moving parts of the slatted frame.
- Check the power cord regularly for damage (abrasions, exposed wires, kinks, pressure marks, etc.) in particular:
  - ▶ after every larger mechanical strain e.g.: crushing with the base frame or rolling over the power cord with an equipment cart.
  - ▶ after relocating/moving and before plugging in the power cable.
  - ▶ during prolonged operating by the user.
- Make sure that the low voltage connection cable of the transformer is not damaged by the base frame or an equipment cart.
- Do not use the area under the slatted frame as a "parking spot" for any utensils.
- Lower the reclining surface in a lower----- position when leaving the patient unattended. This reduces the risk of injury to the patient getting in and out of the bed.
- Keep the hand switch safe from accidentally falling during non-use (Hang it on the hook). Take care that the the cable is not damaged by the moving parts of the slatted frame.

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■ To protect the patient and especially children from unintentionally adjusting the hand switch, place it beyond their reach (e.g. at the foot end) or lock the hand switch when:

- ▶ the patient is not able to control the slatted frame safely or is unable to get out of dangerous positions without help;
- ▶ the patient could be endangered, because of unintentional adjustments to the drive motor;
- ▶ the side rails are fitted:  
(There is a possible danger that limbs could be squeezed when adjusting the back- and thigh rest);
- ▶ unsupervised children are in the room with the slatted frame.

■ Adjustment functions can only be performed by a trained person or in the presence of a trained person.

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■ The reclining surface frame must not dip below the bed box in which it was installed. There is a danger that limbs could be squeezed! When lowering leave 2.5 cm tolerance.

Note the safety instructions in chapters:

- ▶ Safety guidelines for electrical drive systems (see page 26)
- ▶ Cleaning and disinfecting (see page 38)
- ▶ Safety guidelines for exchanging electrical components (see page 55)

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## PRODUCT DESCRIPTION

### NORMAL OPERATING PROCEDURE

■ The slatted frame LIPPE IV, hereafter called slatted frame, serves as a comfortable solution for positioning and easing the care of infirm, frail people in senior or nursing homes. Additionally it was developed as a comfortable solution for in-home care of infirm, handicapped or frail people. It will enhance such this care.

■ The slatted frame is not intended for use in hospitals.

■ This slatted frame can be used for care under the instructions of a physician and can serve in diagnosing, treatment and observation of a patient. It is therefore equipped with a locking mechanism for the electrical adjustment functions.

■ This slatted frame has no special connector for potential equalization. Consider this before combining with additional electrically operated (medical) devices. Further references about additionally applicable precautions can be found:

- ▶ In the instructions for additional electrically operated devices (e.g. compressors for positioning systems, infusion pumps, nutrition probe ...)
- ▶ In the Standard DIN EN 60601-1-1:2002 (Safety of medical electrical systems)
- ▶ In the Standard VDE 0107:1994 (High power systems in hospitals)

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■ The safe working load amounts to 2000 N (~200 kg), whereas a patient weight of 165 kg has been allowed. The remaining 35 kg are allocated for the mattress and attached accessories.

■ This slatted frame is not appropriate for patients less than 150 cm tall. Observe, especially with patients in poor health condition, the safety tips in chapter 4.

■ This slatted frame may only be operated by trained personnel.

■ This slatted frame is useful for multiple reuse. Bear in mind the necessary requirements:

- ▶ See chapter 5;
- ▶ See chapter 6.

**This slatted frame may only be operated under the operating instructions laid out in this instruction manual.**

**Any other use is deemed to be inappropriate.**

**Sale of this product in North America, especially in the United States of America, is prohibited. Distribution and use of the slatted frame, including distribution and use by third parties in these markets is forbidden by the Manufacturer.**

#### SPECIAL FEATURES

- ◆ Electrical height adjustment of the reclining surface from about 29 to 74 cm.
- ◆ Electrical adjustment of the backrest from 0° to about 70°.

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- ◆ Electrical adjustment of the leg rest from 0° to about 40°.
- ◆ Reclining surfaces in various sizes, divided into four areas (see page 5).
- ◆ Trapeze pole with trapeze handle available (accessories).

### STRUCTURAL DESIGN

The slatted frame is shipped disassembled. It can be transported without problems even in apartment homes. It consists of a reclining surface frame in two parts with mounted drive motors for the head and foot sections. The reclining surface frame rests on a base frame. A separate motor controls the up and down movements.

#### Reclining surface frame

The reclining surface frame has four sections: a movable backrest, a firm middle section, a moveable thigh rest and lower leg rest. The back- and thigh rests can be adjusted with the aid of electric motors. The horizontal height of the reclining surface can be adjusted. All adjustment functions are controlled by a hand switch.

#### Base frame

The base frame is equipped with height adjustable pedestal feet leveling on uneven floors. There are four fastening latches within the existing bed box. The base frame can be adapted to different bed widths.

#### Electric adjustment system

The electrical adjustment system of this slatted frame is error protected, flame-retardant (V0) and consists of:

- an external transformer.

The transformer consists of a power cable, the transformer, and a low voltage connection cable.

The transformer creates a 24 volt low voltage, which is safe for patient and user. The

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transformer supplies all drive motors with the 24 volt safety low voltage through cables and a wiring harness. The connector on the base frame is moisture proof.

- a hand switch with a sturdy hook.

The user can lock the adjustment functions of the hand switch if the poor clinical condition of the patient warrants it.

- one drive motor with integrated control for the horizontal height.
- one drive motor with integrated control for the thigh rest.
- one drive motor with integrated control for the backrest.

#### DESCRIPTION OF MATERIALS

The slatted frame is built for the most part from formed steel, the surface of which is covered with a polyester powder coating or has a metallic coating of zinc or chrome.

The spring wood slats consist of wood or wood materials, the surface which is sealed.

**All surfaces are harmless to skin contact.**

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

## ASSEMBLY AND START UP PROCEDURE

The slatted frame is shipped disassembled.  
The assembly is taken care of at the destination by qualified personnel supplied by the operator. Assembly can be done by one or two persons.

The slatted frame ships in two cartons:


**Carton 1**  
Divided reclining surface frame with motor drives, frame connectors, mattress brackets.

**Carton 2**  
Base frame with motor, transformer unit, hand switch and instruction manual.

## ASSEMBLY

### BASE FRAME AND RECLINING SURFACE FRAME

- ☞ Remove the contents of both cartons.
- ☞ Put the base frame [12] in the bed box. The connecting socket [14] must point to the foot end of the bed box.
- ☞ Take the transformer unit [15] out of its wrappings.
- ☞ Pull the plug out of the connecting socket [14].
- ☞ Plug the 6-pole plug (a) of the transformer unit in this connecting socket. The 6-pole plug has a safety groove. Press the plug firmly into the connecting socket.
- ☞ Positioning the transformer unit [15] secure in the bed box.

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☞ Plug the power plug of the transformer into the wall socket.

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The two frame connectors [9] are secured to the top end reclining surface frame before leaving the factory.

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Please note, that the frame connectors for the right and left side are different.

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☞ Remove one frame connector from the reclining surface frame by loosening the knurled screws.

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
**Leave the other frame connector secured in position!**

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☞ Insert the support tube of the removed frame connector into the open end of the V-shaped tube.

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The inside of the V-shaped-tube is lubricated. If you accidentally pull out the frame connector, do not wipe off the lubricant!

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
☞ On the hand switch [10], push on the left button of the middle pair.

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The scissors of the base frame now drive upwards as long as you press the button.

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☞ Move the scissors of the base frame upwards until the V-shaped tube is about 5 cm above the top edge of the bed box.

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
☞ Take the top end of the reclining surface frame [1] and insert the tube support of the frame connector [9] into the open end of the V-shaped tube.

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☞ Replace the opposite frame connector and secure to frame.

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The pin of the frame connector must be located in the hole of the reclining surface frame. Tighten with the knurled screw.

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☞ Fold the thus assembled half of the reclining surface frame back onto the base frame.

☞ Loosen the 4 knurled screws in the top half of the reclining surface frame. Do not unscrew the knurled screws completely, but leave them in the threads for about 2 revolutions.

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☞ Lift the bottom half of the reclining surface frame above the top half of the reclining surface frame.

☞ Insert the halves of the reclining surface frame into each other.

☞ Tighten the 4 knurled screws securely. Do not use pliers to tighten the screws!

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### CONNECTION OF MOTOR DRIVES


☞ Now disconnect the main power plug.

☞ Next connect the drive motors of the head and foot sections with the wiring harness [17].

☞ Insert the plug with the continuous cable into the connector for the motor drive [5] for the head section.

Secure the cable by snapping it tightly into the opening of the strain relief.

☞ Plug the last cable at the end of the wiring harness into the connector for the motor drive [3] of the foot section.

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Secure the cable by snapping it tightly into the opening of the strain relief.



☞ Now plug the power plug of the transformer into the wall socket.

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**Observe without fail:**

When laying the cable in the base frame, make sure that through the up and down movement no shearing or squeezing damages the cables.

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☞ Loosen both knurled screws of the guide rollers. Do not unscrew the knurled screws completely, but leave them in the threads for about 2 revolutions.

☞ Completely insert the two guide rollers into the tube to the limit and leave them there.

☞ Fold the complete reclining surface frame into the horizontal position.

☞ Take the weight of the reclining surface frame and withdraw the guide rollers in turn until the mark appears. Align the guide rollers into the left and right guide bar [18].

☞ Tighten the 2 knurled screws **securely**! Do not use pliers to tighten the screws!

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**Take special note:**

Check the correct position of the guide rollers in the two guide bars. Ensure that both guide rollers are running freely in the guide bars. If necessary adjust sideways equally.

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After the adjustment of the guide rollers firmly **tighten** the knurled screws, otherwise the reclining frame can jump out off the guide bars!

The guide bars are lubricated. Do not wipe off the lubricant!

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### TO PUT THE SLATTED FRAME INTO AN EXISTING BED BOX

☞ Lower the reclining surface frame all the way down in order to ascertain clearance within the bed box. Ensure there is free space equally on all sides.

☞ Raise the reclining surface frame into the highest position.

☞ Use the pedestal feet [13] to adjust to any unevenness of the floor. Turning the feet to the right will raise the frame, turning them to the left will lower it.


The slatted frame can now be secured to the bed box.

☞ Loosen the 4 knurled screws near the fastening latches [11] and withdraw the fastening latches until they make contact with the bed box.

☞ Rotate the fastening latches up and use a wood screw to screw them to the bed box.

Choose the correct length of wood screws.

☞ Tighten the 2 knurled screws near the fastening latches **securely**! Do not use pliers to tighten the screws!

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
The slatted frame is now firmly attached to the bed box.

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In case the bed box requires the base frame to be higher, extenders for the pedestal feet are available.

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Remove the feet from the pedestals, screw in the extenders and replace the feet under them.

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The top of the reclining surface frame [1] must be 60 mm over the bed box!

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The surrounding space between the slatted frame and the bed box must be a minimum of 25 mm to avoid squeezing limbs while reaching into the bed box.

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The mattress brackets [4] should not be below the top of the existing bed box. (Danger of injury!)

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- There must be enough room for the whole assemble of the slatted frame. There must be no furniture, window sills, etc. in the way
- The drive area underneath the lift motor must be kept free.
- A properly installed 230 volt wall plug must be available within the vicinity of the slatted frame.
- Take care when attaching further accessory equipment (e.g compressors for positioning systems, etc) that their attachment is safe and their functionality is not impaired. Especially important is hereby also the safe placement of

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all moveable connecting cables, tubes, etc. If you have any further questions or need clarification contact the manufacturer of the accessory equipment or BURMEIER.

**Observe the following points when setting up the slatted frame to minimize possible fire hazard through outside effects. Instruct the user and the authorized personnel in these points.**

- If possible use only fire-retardant mattresses and bed sheets.
- Only use appropriate firm mattresses according DIN 13014 with a density of 35 to 40 kg/m<sup>3</sup>. Further this mattress should be inflammable according to DIN 597 Section 1 and 2.
- Avoid smoking in bed because of the danger of fire.
- Use only technically faultless appliances (e.g. heating blankets) and other electronic devices (e.g. lamps, radios).
- Take care that these devices are only used for the purpose intended and that they cannot accidentally get on top of or underneath the bed sheets (Fire hazard because of overheating).
- Absolutely avoid the use of extension cords or multi outlet extension cords under the bed (Leaking liquids can cause a fire hazard).
- Always plug the power plug directly into a wall socket. Extension cords and multi outlet extensions should not be used at all.

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## REACHING OPERATIONAL READINESS

■ After successful assembly complete an inspection according to chapter 6.

■ Clean and disinfect the bed before first time use and before every further use according to chapter 5.

**After successfully completing and heeding all steps on pages 17 to 24 the bed is now operationally ready.**

## DISASSEMBLY

☞ Raise the slatted frame to the highest possible horizontal position.

☞ Remove the wood screws from the fastening latches of the bed box.

☞ Loosen the knurled screws on the fastening latches. Push the fastening latches all the way into the base frame. Secure the fastening latches with the knurled screws.

☞ Lower the slatted frame into the lowest horizontal position.

☞ Pull the power plug out of the wall socket.

☞ Lift the slatted frame out of the bed box.

☞ Pull out the plugs of the drive motors for the head and foot sections.

☞ Remove the knurled screws from the frame connectors.

☞ Loosed the knurled screws at the guide bushes.

☞ Lift up the reclining surface frame towards the head end.

☞ Plug the power plug back into the wall socket.

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**SPECIAL SAFETY INSTRUCTIONS FOR THE DRIVE UNIT**

- Take care that any adjustments of the bed no limbs of patient, user, and any other persons, especially children can be jammed and injured by getting under the moveable parts of the slatted frame or in combination with the bed frame.
- When using accessories with the electrically adjustable slatted frame it is necessary to: make sure that the arrangement of these accessories will, when adjusting the back- or thigh rest, not cause any places that might squeeze or hurt the patient. If this is not the case the operator must safely disable the adjustment function of the back- and thigh rests.
  - ▶ To do so place the hand switch out of the patient's reach (e.g. foot section) or lock the adjustment functions on the hand switch.
- Make sure that the power and hand switch cables cannot be jammed or in other ways be damaged.

- Make sure that when adjusting the bed no obstacles such as furniture or slanted ceilings are in the way. This will prevent damage.
  - Make sure that there are no objects on or in the spaces of the base frame.
- This could cause damage to the slatted frame and/or the base frame.
- Take care when laying the hand switch cable that it is not squeezed by any pieces of the slatted frame

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- Make sure that the connector socket is always closed with the plug when it is disconnected from the transformer.

- Take care that the wiring harness or the cable of the hand switch is not pulled or squeezed in any way.

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- ♦ An electronic overload detector switches off the motor drives, when the load is too great. After clearing the overload and pressing the hand switch again the motor drives will work again.

- ♦ A continuous run time of 2 minutes should not be exceeded. Afterwards a pause of no less than 18 minutes should be observed. Another possibility: 1 minute continuous run time and 9 minutes pause.

- ♦ If the maximum continuous run time is not heeded (e.g by constantly "playing" with the hand switch) a thermal safety component switches off the power for safety reasons. Depending on how large the over-load was, it can take several minutes before any more adjustments can be made.

- ♦ The adjustable range of every function is electrically/mechanically only within predetermined ranges possible.

- ♦ Electrical disturbances affecting other nearby electronic appliances (e.g. crackling on the radio) cannot be excluded, despite adherence to all prescribed regulations. In such rare circumstances increase the distance between the appliances, do not use the same power socket or switch the affected appliance off for a while.

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## HAND SWITCH

The electrical adjustments can be made using the hand switch [10].

The adjustable range of every function is electrically/mechanically only within predetermined ranges possible.

Locking mechanisms have been built into the hand switch for safety reasons. For the safety of the patient, adjustment functions can be locked at the hand switch, if, from a clinical point of view and with the approval of the attending physician, the condition of the patient warrants it.

- ◆ The hand switch, with its elastic hook can be hooked onto the bed.
- ◆ The spiral cable allows for greater mobility.
- ◆ The hand switch is water-tight and washable (IP 66).

### OPERATING THE HAND SWITCH

- Train the patient in the functions of the hand switch.
- The motor drives will run as long as the button is held down.
- The adjustment functions, except for the lowered foot position, can be used in both directions.
- Basically the functions of the buttons are as shown :

- ▲ Raise
- ▼ Lower

### Adjusting the backrest

Push this buttons if you want to adjust the backrest angle.

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### Adjusting the reclining surface height

Push this buttons if you want to adjust the height of the reclining surface.

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### Adjusting the thigh rest

Push this button to adjust the thigh rest angle.

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## LOCKING FUNCTION

### Only the operator is authorized to use the locking function!

If the clinical condition of the patient is so critical that any adjustment of the bed by means of hand switch causes danger, the operator must disable the bed immediately. The slatted frame remains in the position it had at the time of disconnecting.

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Do not turn the locking key with force beyond the stop of the lock! The lock or the whole hand switch can be damaged.

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The locking key is delivered attached to the slatted frame with a cable tie on the hand switch. The locking key is not intended for use by the patient. It must be removed from the hand switch. The operator or a person authorized by the physician should keep the locking key.

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■ Turn the respective lock on the hand switch with the locking key clockwise in the locking position.

The color of the respective display will change from green to yellow.

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Lock is vertical

Color of the display is: **green**

Buttons can be pushed ("Clicking sound")

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ca. 15° gedreht

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#### **Drive disabled**

Lock has been turned clockwise by about 15°

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## **LOWER LEG REST**

When the thigh rest is fixed, the lower leg rest **[2]** can be individually adjusted. For this reason there are two adjusting ratchets, which will lock the lower leg rest into position.

#### **Raising**

☞ Raise the lower leg rest on the foot end - not with the mattress brackets - up to the desired level. The lower leg rest will automatically snap into place.

#### **Lowering**

☞ First lift the lower leg rest all the way to the stopping point.

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There is a potential for injury if the lower leg rest falls straight down.

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If you lower the thigh rest, the lower leg rest will be lowered automatically.

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
## EMERGENCY LOWERING OF THE BACKREST

The upright backrest can be lowered manually in case of a power outage or failure of the power supply or the electrical drive systems.

**In this case two users are absolutely necessary!**

**Disregarding these safety guidelines and instructions may cause the backrest to fall out of control, which could lead to serious injuries for the user and the patient!**

Emergency lowering may only be carried out in the case of extreme emergency and by users who have a complete command of the procedures below.

We strongly advise you to practice emergency lowering of the backrest under normal conditions so that in case of an emergency you will be able to react quickly and properly.

- ☞ Before the emergency lowering, ease the load on the backrest.
- ☞ The first user raised the backrest slightly by the outside edge of the head section and hold it in this position.
- ☞ The second user now removes the spring pin from the lift rod. To do this swing away the bent bracket and pull the spring pin along with the bent bracket out of the lift rod of the backrest motor.
- ☞ The motor has now been detached from the backrest and swings down.
- ☞ Now the first user carefully lowers the backrest.

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
**If the backrest falls out of control, the patient and the second user may be injured!**

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☞ To convert the bed back to its normal operating state: Swing up the connecting mount and secure the spring pin, then relocate the bent bracket.

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## MATTRESS BRACKETS

The mattress brackets [7] can be adjusted to accommodate a mattress width of 90 and 100 cm.

☞ Loosen both of the socket head cap screws on the mattress brackets.

☞ Put the mattress brackets onto another pair of holes and tighten the screws firmly.

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## TRAPEZE POLE (ACCESSORIES)


The trapeze pole attached to the bed makes it easier for the patient to get into and out of bed.

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■ The maximum load of the trapeze pole is 75 kg at the front end.  
The trapeze pole is not suitable for rehabilitation training.

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
In both corners of the head section of the reclining surface frame are two trapeze pole sockets [8]. The trapeze pole sockets have a notch at the top (c), which together with the metal pin (b) restricts the swivel range of the trapeze pole. The trapeze pole should be installed on the side where the patient gets in and out of bed.

#### Attaching

☞ Insert the trapeze pole into the socket. The metal pin (b) must be in the notch (c) of the socket.

#### Removal

☞ Pull the trapeze pole straight up out of the socket.

#### Trapeze pole swivel range

Do not swivel the trapeze pole outside of the bed. There is a danger that through the pull, the bed will tip over.

The metal pin must always be in the notch of the socket.

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
## TRAPEZE

A trapeze (grab handle) can be attached to the trapeze pole, it is enclosed with the trapeze pole. The patient can hold on to the trapeze to make it easier to change to a different position.

Regularly check the trapeze at the strap for any damage (see Chapter 6). A damaged trapeze or strap must immediately be exchanged.

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### Durability

There is a date clock on the trapeze. With normal use the trapeze will last for at least five years. Afterwards a sight check is to be performed every six months to determine if further use is permissible.

### Adjustable range of the trapeze

The height of the trapeze can be adjusted by means of the strap to accommodate a range of about 55 cm to 70 cm (measuring from the top of the mattress).

- ☞ Push the loop of the trapeze [14] over the first stud of the trapeze pole.
- ☞ Test by tightly pulling the trapeze downward that it is secure.

**Useful Tip:** The maximum lifting load of the trapeze pole is 75 kg at the front end.

- The height of the trapeze can be adjusted by means of the strap.
- Make sure the the strap is fed correctly through the buckle.
- Make sure that the end of the strap is at least 3 cm above the buckle.

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### MATTRESS BRACKETS LONG (ACCESSORIES)

If you are using a very high mattress, the mattress brackets delivered can be exchanged for extra long mattress brackets.

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## SIDE RAILS (ACCESSORIES)

Side Rails conveniently prevent the patient from falling out of bed accidentally. They do not, however, prevent a patient from deliberate leaving the bed.

### Special safety guidelines for side rails

- Only use technically perfect, undamaged side rails with certified gap dimensions, which can safely be attached.
  - There is a danger that between the side rails and the bed frame, the long opening can cause squeezing.
  - Make certain that only suitable Original BURMEIER side rails are used, which have been approved for use with the bed model you are using or that are already integrated in the bed.
  - Before using side rail beams evaluate and take into account the clinical condition and any characteristics of the build of the patient.
    - ▶ If, for example, the patient is very bewildered or very anxious, forgo the side rails and use alternative safety measures such as a restraining belt or restraint sheets.
    - ▶ For especially small, frail patients, additional protective measures for reducing the distance between bars on the side guard may be necessary. In these cases use side guard foam covers, retraining belts, etc.
- This is the only way to effectively ensure patient safety by reducing the risk of becoming trapped or slipping through.
- Only use appropriate, firm mattresses accord to DIN 13014 with a density of 35 to 40 kg/m<sup>3</sup> and a height of 10 to a maximum of 12 cm. Further this mattress should be inflammable according to DIN 597 Section 1 and 2.

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- If elevated patient positioning systems are employed, such as anti-decubitus mattresses, an effective side guard height of at least 22 cm above the unweighted mattress must be assured. If this requirement is not met, the operator must take suitable measures to safeguard the patient.

- To safeguard the patient against unintentional adjustments, always lock the electrical adjustment of the back- and thigh rests on electrical slatted frame at the hand switch, when the side guards are raised.

Otherwise there is a risk of crushing or trapping limbs when the back- and thigh rest are adjusted.

**If these safety warnings are not observed there is a danger for the patient of being injured by the side rail, of falling out of bed because of too large a distance between the bars or because of damage, the use of unsuitable side rails or the incomplete locking of the side rails.**

- Remove the mattress brackets [7] on the side where the side rail will be installed.

- Attach the side rails in such a way on the reclining surface frame that on the head and foot areas there is a minimum of 25 cm space to the bed frame.

- The side rails click audibly into place in the reclining surface frame.

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- Test the secure and firm placement of the side rails.
- Before lowering the side rail beams make sure that there are no limbs of the patient in the area of the side rail beams. This way you can avoid crushing or injuring the patient.
- When attaching the side rails the reclining surface area extends an additional 3 cm. Check to see if the side rail when lowered sets on the bed box!

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## CLEANING AND DISINFECTING

### GENERAL

Cleaning is the most important requirement for a successful chemical disinfection.

Generally a routine cleaning of the bed, if used by the same patient, is hygienically sufficient

Disinfection of the bed frame is only necessary in case of visible contamination with infectious or potentially infectious material (blood, stool, pus) or by presence of an infectious disease under the direction of a physician.

**The slatted frame must be cleaned and washed down with disinfectant in case of a patient change!**

#### Before beginning the cleaning, please note:

- Unplug the power cable from the wall socket and store it in such a way that it does not come into excessive contact with water or other cleaning solutions (store in plastic bag).
- Make sure all plugs at the drive motors are properly connected.
- None of the electrical components must show any visible damage; otherwise water or cleaning solution could penetrate. This could lead to malfunction or damage to the electrical components.
- Before operating the bed again, make sure that there is no residual moisture on the electrical contacts, by drying or blowing on the electrical plugs.

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- The electrical components must not be cleaned with a water jet, a high pressure cleaner or other similar device! Clean only with a moist towel.
- If you suspect that water or any other form of moisture has penetrated the electrical components, pull the power plug out of the wall socket immediately. If already disconnected from the power supply, make sure it is not plugged in again. Clearly mark the slatted frame as "Defective" and take it immediately out of service. Report this to the operator promptly.
- Failure to follow these directions could lead to major damage of the slatted frame and its electrical components and lead to subsequent malfunctions.

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## CLEANING AND DISINFECTION PLAN

- Remove the bed sheets and give them to the laundry.
- Clean all surfaces including the slatted boards and the reclining surface made of plastic inserts or wire mesh base with a mild and environmentally safe cleaning solution. This also applies to the hand switch.
- If the bed has visible contamination e.g. infectious or potentially infectious material have polluted if, it must be disinfected by wiping after cleaning. Use approved disinfectants for each surface according to the products listed by the DGHM (German Organization for Hygiene and Microbiology). This also applies to all beds from patients with illnesses subject to report according to §6 of the infection control law (*German*

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*abbreviation:* IfSG) Protection against Infection Act, colonization or infection with treatment resistant pathogens (e.g. MRSA, VRE) and all beds from the intensive care and infection control units. It is necessary to use the DGHM approved concentration.

■ Disinfecting the pedestals is only necessary when visibly contaminated by infectious or potentially infectious material.

**Notice:** Continuous disinfection is only necessary for patients with a treatment resistant pathogen (e.g. MRSA) within the hospital.

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## INSTRUCTING THE USER AND QUALIFIED PERSONNEL

To guarantee the correct procedure for cleaning and disinfection is followed, we recommend to correctly instruct the handler and qualified personnel. Thereby communicating that the following points need to be observed:

■ The clean slatted frame is to be transported to the patient's home in such a way that it does not become polluted or contaminated.

■ When disassembling the slatted frame we recommend cleaning and disinfecting it right away. The qualified personnel should be familiar with the steps necessary to thoroughly clean and disinfect (the operator will provide the individual steps necessary). It is necessary to use only the DGHM approved disinfectant in the approved concentration. The disinfectant must be appropriate for the surface being cleaned.

■ The qualified personnel must be equipped with water-tight one-time-use aprons and gloves.

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■ Only fresh and clean towels must be used to clean and disinfect, these must then be given to the laundry.

■ After cleaning the slatted frame, the qualified personnel must also disinfect their hands before going to the next patient. An appropriate disinfectant dispenser should be part of the qualified personnel's gear.

The immediate cleaning of the slatted frame while disassembling has the benefit of not transporting "dirty" slatted frames or frame parts along with clean ones in the same truck. A transfer of potentially infectious spores which may be clinging to the used slatted frame may be avoided. A transfer of spores of a hospital-acquired infection will be avoided with strict adherence to these recommendations.

If the slatted frames are to be stored for a while, they should be covered to keep them free of dust, accidental pollution and contamination.

## CLEANING SUPPLIES AND DISINFECTANTS

Make use of these recommendations to keep the slatted frame serviceable for a long time.

■ The surfaces must be intact; any damage should be repaired immediately.

■ We recommend a (damp) wash down. When choosing the cleaning solution make sure

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reparation is mild (gentle to skin and surfaces) and environmentally friendly. In general a regular household cleaner can be used.

■ For cleaning the coated metal parts use a wet cloth with regular household cleaners

■ The disinfectant used for the wash down should be one approved for the purpose and listed by the DGHM (in the recommended dosage). Disinfectants and concentration listed by the Robert-Koch-Institute (the former Bundesgesundheitsamtes) should only be used in case of an epidemic according to the infection safety law and only on the orders of the appropriate health department.

■ Despite the excellent mechanical resistance of the bed, scratches, markings, etc. which penetrate the coating should be resealed using a suitable medium to prevent the penetration of moisture. Contact BURMEIER or your specialized dealer with any questions.

■ Do not use scouring agents, stainless steel cleaners, abrasive cleaning agents or scouring pads. These substances can damage the surface.

#### Notice:

Surface disinfectants based on aldehyde generally have the advantage of working on a large array of pathogens, have a relatively low protein error and are environmentally friendly. Its main disadvantage is the allergy-causing and irritation potential. Compounds on the basis of gluco-

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protamine do not have this disadvantage and are very effective, though mostly more expensive. Disinfectants on the basis of chlorine splitting combinations, with prolonged contact or in too high doses, can be corrosive to metals, plastics, rubber and other materials. In addition to this, the compounds have a high protein error, are irritants to the mucous membranes and are not environmentally friendly.

#### CONTACT WITH DISINFECTANTS

- Follow the guidelines for each product and manufacturer. Observe the exact dosage! We recommend automatic dosage apparatus.
- Always mix the solution with cold water, to avoid vapors that aggravate the mucous membranes.
- Wear gloves to avoid direct skin contact.
- Do not store the mixed disinfectant in open containers with the cleaning rag inside. Close the container!
- Use closeable bottles with dosage pumps for wetting the cleaning rag.
- After finishing the surface disinfection ventilate the room.
- Use a wipe down disinfecting method, do not use a spray! Using the disinfectant as a spray causes the disinfectant to turn into a mist which could be inhaled.
- The wipe-down disinfection plays a major part.
- Do not use alcohol for larger areas.

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Operators of nursing beds are obligated to keep medical products in a safe state throughout their entire service life. This also includes regularly carrying out expert maintenance as well as safety checks.

### Notice for Operators

This slatted frame has been constructed and built to last safely for a long time. With appropriate operation and use, this slatted frame has a life expectancy of 2 to 8 years. The life expectancy is based on the conditions and frequency of use.

Through repeated transport, assembly and disassembly, improper operation as well as long time use, it is possible that damage, defects and deterioration can occur. These deficiencies can lead to danger if they are not recognized and quickly fixed.

Because of this there are legal bases for the regular enforcement of tests to ensure the safe condition of this medical product.

According to § 4 of the Medizinprodukte-Betreiberverordnung operators are obliged to maintain the product. Because of this the following regular inspection and function controls are to be

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performed by the operator as well as the user.

**Point out these necessary tests to the user.**

## BY THE USER

Along with the regular substantial tests by technically qualified personnel, the normal user (nurses, supporting relatives, etc.) must also perform a sight and function check at shorter, regular intervals, as well as before a new patient.

- If damage or a malfunction is suspected, take the slatted frame out of service immediately and disconnect it from the power until an exchange or repair can be made!

- Contact your operator regarding exchange or repair of faulty parts.

### Recommendation:

All electrical and mechanical components should be tested monthly. The main power cable as well as the hand switch cable should be checked after every mechanical load and after every relocation.

For this purpose use the checklist on the following **page 47**.

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**Checklist: Test protocol by the user**

Test	OK	Not OK	Defect description
<b>Visual check of electrical components</b>			
Hand switch and cable	Damage, Cable Routing		
Power cord	Damage, Cable Routing		
Hand switch	Damage, Foil		
Transformer unit	Damage, Cable Routing		
<b>Visual check of mechanical components</b>			
Trapeze pole/Socket	Damage, Deformation		
Bed frame	Damage, Deformation		
Spring wood slates	Damage, Splinters		
Bed box	Damage, Splinters		
Reclining surface frame	Damage, Deformation		
Side rail	Damage, Splinters		
<b>Function test of electrical components</b>			
Hand switch	Functiontest, Blocking function		
Strain relief, Power cable of transformer unit	Tight fit		
<b>Function test of mechanical components</b>			
Emergency lowering of backrest	Test according to instruction manual		
Socket head cap screw	Tight fit		
Knurled screws	Tight fit		
Side rail	Safe locking/unlocking		
Lower leg rest	Latching		
Accessories (e.g. trapeze pole, trapeze)	Secure, Damage		
Signature of the tester:	Result of the test:		Date:

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## BY THE OPERATOR

The operator of the slatted frame according to MPBetreibV § 4 is obligated with every new assembly, every maintenance and throughout the operation, to perform regular inspections, to ensure the safe condition of the slatted frame.

These tests are to be repeated during regular maintenance and depending on the use according to MPBetreibV § 4 and from the Directive of the German Employers Liability Insurance Association according to BGV A2 (former VBG4).

- ◆ Observe this sequence for the inspection according to VDE 0751:
  - I. Visual test
  - II. Electric Measurements
  - III. Function Test
- ◆ The performance of the function test as well as the rating and documentation of the test results can, according to MPBetreibV § 4 only be performed by experts with the appropriate knowledge, training and experience as well as the tools necessary to perform it according to the rules.
- ◆ The performance of electrical measurements according to VDE 0751 can, with availability of proper measuring equipment, also be performed by an electro-technically trained person (according to BGV A2 (formerly VBG 4)) with medical and appliance specific training.
- ◆ The assessment and documentation of the test results can only be performed by an expert with the corresponding knowledge, training and experience.
- ◆ Testing of the transformer requires only measurement of the output.  
Thus, after replacing the transformer with a new unit, the bed can be put to use.

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
♦ BURMEIER offers transformer testing as customer service. Send the transformer unit back to BURMEIER and receive a tested replacement unit. Contact us for further details; Address see page 55.

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#### Preperation:

- Unplug the power plug of the transformer from the wall socket.
- Unplug the 24 volt cable from the connecting socket.
- Plug the plug of the 24 volt cable in the adapter (extra accessory; request from BURMEIER).
- Connect the adapter with the test probe of the test unit.
- Plug the power plug of the transformer into the power socket of the test unit.

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### Test method:

- Output power test: direct or differential according to DIN VDE 0751-1:2001
- Output power test according to the test unit manufacturers instructions.

### Limit:


- Output power less than 0.1mA.

### Test cycle:

We recommend that the yearly inspection, which, depending on usage, with accountability according to demonstrably observing the 2%-error-quota (see BGV A3: § 5, table 1 B) can maximally be extended to a total of two years.

If damage or a function error is suspected take the slatted frame out of service immediately and disconnect it from the power until an exchange or repair can be made!

The test protocol draft on **pages 53 and 54** can be used.

Model	Bettenmodell
Type	Artikelnummer
Serie	Seriennummer

Model	Bed modell
Type	Item Number
Serial	Serial Number

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## SPARE PARTS

The corresponding replacement parts can be obtained from BURMEIER by specifying the item number, the order number and the serial number. The necessary details are found on the type plate under the cross tubing of the reclining surface frame.

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<b>Customer / med. Facility / practice:</b>					
<b>Address:</b>					
Carried out: <input type="checkbox"/> Repeat inspection		<input type="checkbox"/> Inspection prior to initial operation			
<input type="checkbox"/> _____		<input type="checkbox"/> Inspection following repairs / servicing			
<b>Equipment type:</b> <input type="checkbox"/> Hospital bed <input type="checkbox"/> Care bed <input checked="" type="checkbox"/> Slatted frame		<b>Protection classification:</b> <input type="checkbox"/> I <input checked="" type="checkbox"/> II			
<b>Bed type:</b> LIPPE VI		<b>Inventory Number:</b>			
<b>Location:</b>		<b>Serial Number:</b>			
<b>Manufacturer:</b> Burmeier GmbH & Co. KG, Germany		<b>User-specific parts:</b> none			
<b>Testing equipment used (Type / Inventory Number):</b>		1.			
<b>MPG Classification:</b> Class I		2.			
<b>I. Visual Inspection</b>			<b>OK</b>	<b>Not OK</b>	<b>Defect Description</b>
<b>Visual inspection of the electrical components</b>					
Stickers und type plates	Presence, legible				
External transformer housing	Damage				
Motor housing and lifting rods of motors	Correct position, Damage				
Hand switch	Damage, Foil intact				
Motor-, hand switch- and power cables	Damage, Routing				
<b>Visual inspection of the mechanical components</b>					
Stickers and type plates	Presence, legible				
Trapeze pole/socket	Damage, Deformation				
Bed frame	Damage, Deformation				
Spring wood slates	Damage, Splinters				
Reclining surface	Damage, Deformation				
Wooden boards	Damage, Splinters				
Welded seams	Split welded seams				
Side rail	Damage, Splinters				
Socket head cap screw	Tight fit				
Knurled screws	Tight fit				
Wearing parts such as joints	Damage				
<b>II. Electrical measurement of transformer according to VDE 0751-1:2001-10 (Output power test: direct or differential according)</b>					
Proceed as follow:			<b>OK</b>	<b>Not OK</b>	<b>Defect Description</b>
1. Plug the transformer power cable into the test socket on the measuring instrument.					
2. Insert the plug of the 24 volt cable into the special adapter (extra accessory; request from BURMEIER).					
3. Connect the special adapter with the test probe of the test unit.					
4. Start procedures with the measuring instrument.					
	Limit Value	Value of first measurement	Current actual value		
Result: Bed SK II (Type B)	0,1 mA	mA	mA		

III. Performance check		OK	Not OK	Defect Description
<b>Performance check of the electrical components</b>				
Motion limit switching of the motors	Automatic cut-off			
Hand switch, blocking functions	Test according to instructions for use			
Drive motors	Noise level			
Drive motors	Test according to instructions for use			
Strain relief and security of routing the plug	Secure			
Strain relief for the transformer cable	Secure			
<b>Performance check of the mechanical components</b>				
Joints and pivots	Smooth operation			
Emergency lowering of the backrest	Test according to instructions for use			
Side rail	Safe locking/unlocking			
Lower leg rest	Latching			
Accessories (e.g. trapeze pole, trapeze)	Secure, Damage, Suitability			
<b>Result of the inspection:</b>				
<b>All values within permissible range:</b>		<input type="checkbox"/> yes <input type="checkbox"/> no		<b>Next inspection date:</b>
<b>Inspection passed:</b>		<input type="checkbox"/> yes <input type="checkbox"/> no		
<b>If inspection was not passed:</b>		<input type="checkbox"/> Defective, do not use the bed ! => Repair <input type="checkbox"/> Defective, do not use the bed ! => Take out of Service <input type="checkbox"/> Bed does not meet the safety standards		
<b>Test approval sticker applied:</b>		<input type="checkbox"/> yes <input type="checkbox"/> no		
<b>Part of this inspection report:</b>				
<b>Comments:</b>				
<b>Inspected on:</b>	<b>Inspector:</b>		<b>Signature:</b>	
<b>Evaluated on:</b>	<b>Operator/Expert:</b>		<b>Signature:</b>	

Test protocol of  
an inspection of  
electro-medical  
equipment  
according to  
DIN VDE 0751-1:  
2001-10  
Page 2 of 2

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## MANUFACTURER'S ADDRESS

In order to maintain operational safety and the right to claim under warranty, only original BURMEIER replacement parts may be used!

For ordering replacement parts, customer service requests and other questions, please contact:

**Burmeier GmbH & Co. KG**  
(A subsidiary of the Stieglmeyer Group)

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## REPLACEMENT OF ELECTRICAL COMPONENTS

### Dangerous voltage, electric shock possible

- Before doing any work unplug the main power cable from the wall socket!
- Any work on the electrical equipment may only be performed by customer service, the motor manufacturer or qualified and authorized electricians under consideration of the applicable VDE rules and safety guidelines!

- Only dismantle the control unit and electrical motors when the slatted frame is in the horizontal position. Otherwise there is the risk of being injured by falling parts of the mattress base.

- All the drive components are maintenance-free and must not be opened. In the event of a malfunction, the corresponding components should always be replaced in full!

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### TROUBLESHOOTING

The following table on page 58 is a guide to rectifying common malfunctions:

If a function fault is detected during operation, which the user cannot rectify with the use of the trouble-shooting table aid, the specialized staff for maintenance and repair of the appropriate owner is to be called.

The user may **under no circumstances** try to eliminate defects at the electrical equipment!

#### Dangerous voltage, electric shock possible!

Any work on the electrical equipment may only be performed by service personnel, the motor manufacturer or qualified and authorized electricians as required by the applicable VDE rules and safety guidelines!

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### Trouble-shooting table aid

Problem	Possible causes	Solution
Hand switch / drive system not functioning	<ul style="list-style-type: none"> <li>■ Hand switch plug not properly connected</li> <li>■ Transformer plug (24 volt) not properly connected</li> <li>■ Hand switch or drive system defective</li> <li>■ Functions are blocked on hand switch</li> </ul>	<ul style="list-style-type: none"> <li>☞ Test the plug connection</li> <li>☞ Plug the transformer plug into the connection socket</li> <li>☞ Inform your operator about necessary repairs</li> <li>☞ Unlock functions (see page 29)</li> </ul>
Green LED on transformer not lit and drive systems are without function	<ul style="list-style-type: none"> <li>■ Power cable not properly connected</li> <li>■ Wall socket without power</li> <li>■ Transformer defective</li> </ul>	<ul style="list-style-type: none"> <li>☞ Plug in the power cable properly</li> <li>☞ Test the wall socket and fuse box</li> <li>☞ Exchange the transformer; inform your operator about necessary repairs</li> </ul>
Green LED on transformer lit continuously, yellow LED not lit when hand switch button is pushed	<ul style="list-style-type: none"> <li>■ Hand switch is defective</li> <li>■ Transformer plug (24 volt) not properly connected</li> <li>■ Functions are locked on hand switch</li> </ul>	<ul style="list-style-type: none"> <li>☞ Exchange the hand switch</li> <li>☞ Test the plug connection</li> <li>☞ Unlock the functions (see page 29)</li> </ul>
Drives operate in one direction only	<ul style="list-style-type: none"> <li>■ Hand switch or drive system is defective</li> </ul>	<ul style="list-style-type: none"> <li>☞ Inform your operator about necessary repairs</li> </ul>
Drives stop suddenly after prolonged adjustment times	<ul style="list-style-type: none"> <li>■ Heat switch in transformer has been activated</li> </ul>	<ul style="list-style-type: none"> <li>☞ Pay attention minimum pause (see page 27)</li> <li>☞ Exchange the transformer; inform your operator about necessary repairs</li> </ul>

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**Benennung**

## ACCESSORIES

Only the Original BURMEIER accessories listed below may be used. If other accessories are used, we can not be liable for any resulting injuries, defects and dangers.

**Description**

Trapeze pole with trapeze
Side rail bars
Mattress brackets, long (includes mounting hardware)

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When using side rail bars, drip/I.V., etc. with electrically adjustable beds, please note:

Make sure that the arrangement of accessories does not produce any crush or shearing zones for the patient when the back- and leg rests are adjusted. If this cannot be ensured, the user must safely prevent the patient from adjusting the back- and leg rest. For this, lock the hand switch.

► To do so, attach the hand switch out of reach ( e.g. at the foot end of the bed), or lock hand switch adjustments.

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**Assembled slatted frame:**

Total weight (size 200 x 90 cm)
Safe load
Patient weight, max.

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**Disassembled bed:**

Reclining surface frame
Base frame
Trapeze pole with trapeze handle

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Operating time
----------------

Protection, splash proof
--------------------------

Classification, Protection class II, Type B, MPG classification Class I, not for use in explosive atmospheres.

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	<b>Motor for reclining surface height adjustment</b>
	Type
	Power / Size / Lift
	Input Voltage
	Operating Time
	Protection

#### **Motor for backrest**

	Type
	Power / Size / Lift
	Input Voltage
	Operating Time
	Protection

#### **Motor for thigh rest**

	Type
	Power / Size / Lift
	Input Voltage
	Operating Time
	Protection

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**IP 44**

#### **Explanation of Icons :**

♦ Unit with Type B application according to IEC 601-1 (Special protection against electric shock)

♦ Protection Class II device, double insulation

♦ Transformer with thermal fuse block

♦ Safety transformer in accordance with VDE 05510

♦ Caution! Consult manual

♦ The electrical equipment offers protection from internal dust deposits and is splash-proof from all sides

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♦ The electrical equipment is dust-proof from internal dust deposits and is splash-proof from all sides

♦ Conformity mark in accordance with Medical Products directive 93/42 EWG

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## AMBIENT CONDITIONS

### Noise level during motor adjustments

The ambient conditions stated below must be observed:

#### In storage

Storage temperature  
Relative humidity

#### In operation

Ambient temperature  
Relative humidity  
Air pressure

## APPLIED STANDARDS AND GUIDELINES

Medical Products Directive
Repeat testing ... of medical electrical units
Technical aids for handicapped persons
Handicapped friendly construction
Risk analysis for medical products
Adjustable beds for handicapped persons
Safety for electro-medical equipment
Electromagnetic tolerances
Electrically operated hospital beds (adapted for application for in home care)

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Classified as class I active medical product( in accordance with MPG § 13).

Fulfills all additional safety requirements by the Highest State Authority from May 22, 2001 for nursing beds.

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### DISPOSAL INSTRUCTIONS

- ◆ This slatted frame is - if electrically adjustable - classified as a commercially used electric appliance (b2b) in accordance with WEEE-Guidelines 2002/96/EG (Electronics Law).
- ◆ Exchanged electronic components (drives, control unit, hand switch, etc.) of these slatted frames are considered electronic scrap and must be disposed of according to the WEEE-Guidelines.
- ◆ With beds that came into use after August 13, 2005 (see designation "Serial" on the type label at the head of the bed), the operator is required to not dispose of the electrical components in the municipal waste area, but to send them directly to the manufacturer. BURMEIER and their service and sales partners will takes these parts back.
- ◆ Our general terms and conditions apply for this redemption.
- ◆ No longer useable rechargeable batteries have to be disposed of according to regulations and do not belong in the house garbage.
- ◆ With all disposed components the owner must make sure that they are not infectious or contaminated.
- ◆ In the event that the slatted frame must be scrapped, the plastic and metal parts used must be disposed of separately, following the proper disposal procedure.

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hereby declare that the product named below:

### Slatted frame LIPPE IV

complies with the regulations of EC Guidelines 93/42/EWG for medical products and fulfills all additional safety requirements by the Highest State Authority from May 22, 2001 for care beds.

Applied Standards:

- |                              |   |
|------------------------------|---|
| ● EG-Richtlinie 93/42        | Guideline for medical products  |
| ● VDE 0751-1:2001-10         | Repeat testing of ... medical electrical units                                    |
| ● EN 12182:1999-11           | Technical aids for handicapped people   |
| ● DIN 32977-1:1992-07        | Handicapped fair constructions  |
| ● DIN EN 1970 : 2005-10      | Adjustable beds for handicapped people<br>(as far as product specific)            |
| ● DIN EN 60601-1: 1996-03    | Safety for electro-medical equipment  |
| ● EN 60601-1-2:2002-10       | Electromagnetic tolerances  |
| ● DIN EN 60601-2-38/A1: 2001 | Electrically operated hospital beds<br>(reduced to applications for in home care) |

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No part of this manual may be reproduced with-  
out the prior written permission of the publisher!

All rights reserved.  
Subject to technical changes!

178691  
Doc-Item: Lippe4\_v2  
Updated: February 26, 2007

Realisation: [info@docnik.de](mailto:info@docnik.de)